

Nailexpert by **WARTNER**[®]

Medical device for the elimination of nail fungus

Use of Nailexpert by Wartner

Nailexpert by Wartner is a medical device for the elimination and prevention of nail fungus or nails with inclination to fungal infections (both yeasts and moulds). It is designed for treatment of both toe- and fingernails. The intended use is for the treatment of light to medium topical onychomycosis. This preparation is highly recommended to patients who are often exposed to a mycosis environment (frequent visiting of swimming pools, sport places, saunas, public showers etc.). Nailexpert by Wartner is intended for adults only.

Treatment method

Nailexpert by Wartner is a medical device offering a new targeted treatment using a unique formula and application brush. Nailexpert is able to eliminate nail fungus fast through:

1. Preventing growth & further spread (with the user or other persons) of the infection

- by creating an advanced film which avoids infection of other nails and forms an inhospitable environment for the fungus

2. Facilitating regeneration towards a healthy nail

- by the moisturizing & nourishing features of its active ingredients

The mode of action to prevent further spread and growth of infection is based on film-forming silicone oils. The selected silicones have film forming properties and create a non-permeable film on nails. This film defends nails against onychomycosis and other dermatophytes. It also inhibits already present fungi from accessing to oxygen and other nutrients, essential for their growth and life (extracted from living skin and nails) by covering their surface by lipophilic and hydrophobic film.

In addition, to facilitate regeneration of the nail, Nailexpert by Wartner contains D-panthenol, which has regenerating and, in combination with glycerine, also moisturizing & nourishing properties to strengthen the regeneration towards a healthy nail.

Product application

- 1) Prior to application wash your hands or feet carefully, preferably with soap-free mild surfactant preparation. Before use shake properly.
- 2) Preparation should be applied on clean, nail-polish-less, pared nails (without any varnish remover residue).
- 3) Apply a sufficient amount of the preparation with the foam applicator to the affected area. This sufficient dosage is ensured by the applicator.
- 4) Apply to the affected area and let it dry. Repeat 2 times.
- 5) Apply twice a day
- 6) During the treatment period do not use any other preparation as nail polishes, powders, sprays, etc. otherwise efficacy of preparation might be reduced.
- 7) Treatment should last until the total nail fungus symptoms are reduced or until the packaging is used up. Treatment period depends on the nail growth, however the usual treatment period lasts 4-12 weeks. First visible results after 4 weeks.

400 applications included.

Precautions

Preparation is designed for the treatment of adults

Do not use on open wounds

Do not use in case of damaged packaging

For intended use only

For use in other cases than indicated (e.g. for children), a medical professional should be consulted

In case of accidental ingestion of preparation by children, a medical professional should be immediately consulted and the packaging of the preparation should be shown

Do not use after expiry date indicated on the packaging!

Keep out of reach of children

In case of any uncertainty, a medical professional should be consulted and the packaging of the preparation should be shown

Contraindication

The product is generally very well tolerated. Interaction with other preparations or drugs has not been observed. Nevertheless, it is recommended not to combine with other preparations with same indication unless consulted with a medical professional. In case of any side effects or any abnormal reactions, a medical professional should be immediately consulted and the packaging of the preparation should be shown.

Ingredients: Aqua, Cyclomethicone, Glycerin, D-Panthenol, Sodium Dihydroxycetyl Phosphate, Sodium Bicarbonate, Dimethiconol, Alcohol denat., Lactic Acid, Hydrogen Peroxide, C.I. 77 820, C.I. 77 400

Storage

Keep in a dark and dry place at room temperature.

Content

- Tube with applicator 4 ml
- Instructions of use leaflet

Clinical tests of preparation proved the following after 4 weeks of use:

- Nail surfaces are visibly improved (smooth, reduction of dryness, friableness and fragility)
 - Nails are softer and can be adjusted easier filing, cutting.
 - Relief from itching
- Patients involved in testing applied the preparation strictly according to the recommended indications. Product consumption varied depending on the size of treated area.

For more information: www.wartner.com



Omega Technika Farnham Drive, Finglas, Dublin 11, Ireland
Distributor: OMEGA ALTERMED s.r.o., Tomáškova 30,
821 01 Bratislava, SK, www.omega-altermed.eu

4 ml CE₁₆₃₉



Latest revision date:
February 2011